

Alaris pump
primary - drug library

Step	Failure Modes	Severity	Probability	Score	Detectable	Prevent?	CAUSES						ACTIONS							
							Training	Policy & Proc	Technology	Environ	People	Training	Policy & Proc	Technology	Environ	People				
1	Turn system on	not done	L	L	L	Y	N													
		dead battery	L	L	L	Y	N													
		broken button	L	L	L	Y	N													
2	System self-tests	total breakdown	L	L	L	Y	N													
3	Select yes/no for new pt	lack of knowledge of pt status	M	L	LM	N	Y	X												
									X											
										X										
4	Select profile	choose incorrect profile	M	H	MH	Y	Y	X	X											
									X	X										
									X	X										
		repetitive button pushes to find correct profile	L	H	M	Y	N													
		can't select new profile while pump is running	H	H	M	Y	Y				X									
5	Select channel	incorrect button (brain vs muscle)																		
		tubing crossed	H	H	H	N	Y	X	X											
									X	X										
										X										
		chart channel label, rather than infusion name (labels ABC right to left, even if channel added)								X										
6	Select drug library	select basic infusion instead																		
7 I	Select drug	Incorrect drug	H	H	H	Y	Y			X										
		incorrect conc	H	H	H	Y	Y			X										
								X	X											
	Drug calc	enter manually calculated rate	H	H	H	N	Y			X										
										X										
											X									
												X								
		incorrect weight	M	H	MH	N	Y			X										
										X										
		vol entered is not vol in bag	H	H	H	N	Y	X	X											
7 II	Select dose calc	no drug label	M							X										
		incorrect weight	M	H	MH	N	Y													
		enter rate in volume field																		
		enter volume in rate field																		
		double key stroke - skipped screen	M	H	MH	N	Y	X												
8	Enter volume to be infused	incorrect key pads	M	H	MH	Y	Y													
		misread order	M	H	MH	N	Y				X									
										X										
		vol entered is not vol in bag	H	H	H	N	Y	X	X											
9	Review verification screen & correct if needed	OK without reading data on screen	M	M	M	N	Y	X	X											
		read incorrectly	M	M	M	Y	Y													
10	push start																			